

Reproducible bioinformatics at scale: **nf-core + Nextflow**

Phil Ewels · Senior Product Manager for Open Source, Seqera
Bioinformatics CRO Webinar · February 2026

phil.ewels@seqera.io
 ewels.bsky.social  ewels

ABOUT ME

Phil Ewels

Senior Product Manager, Seqera

Creator of MultiQC. Former head of genomics development at Swedish National Genomics Infrastructure - SciLifeLab, Sweden. Co-founder of nf-core and PM + contributor to Nextflow.



10+

Years in bioinformatics



nf-core

Co-founder, core team member



MultiQC

Creator & maintainer



THE PROBLEM

The reproducibility crisis

70%

of researchers failed to reproduce
another scientist's experiments

Nature, 2016 · 10.1038/533452a

50%

failed to reproduce
their own experiments



Complex
dependencies



Diverse compute
environments



No workflow
standards



Poor data
provenance



Workflow orchestration for reproducible scientific pipelines

Flow-based syntax

Simple, readable syntax for complex data pipelines

Resume & cache

Automatic caching — resume from any failure point

Software built-in

Docker, Singularity, Conda for full reproducibility

Scalable

One sample to thousands with auto-parallelization

Multi-platform

Laptop, HPC cluster, or any major cloud provider

Active ecosystem

140+ curated pipelines via the nf-core community



HOW IT WORKS

Four simple steps

1.

Define

Write processes (commands) in Nextflow code

2.

Connect

Link processes via channels. Nextflow handles data flow automatically

3.

Configure

Set compute environment, containers, and resources

4.

Execute

Run on laptop, HPC, or cloud

```
process {  
  
    input:  
    path 'sample.fq.gz'  
  
    output:  
    path 'report.html'  
  
    script:  
    """  
    mytool sample.fq.gz \  
        -o report.html  
    """  
}
```



HOW IT WORKS

Four simple steps

1. Define

Write processes (commands) in Nextflow code

2. Connect


Link processes via channels. Nextflow handles data flow automatically

3. Configure

Set compute environment, containers, and resources

4. Execute

Run on laptop, HPC, or cloud



HOW IT WORKS

Four simple steps

1. Define

Write processes (commands) in Nextflow code

2. Connect


Link processes via channels. Nextflow handles data flow automatically

3. Configure

Set compute environment, containers, and resources

4. Execute

Run on laptop, HPC, or cloud



KEY CONCEPTS

Nextflow key points

Language

Code for describing workflow steps and data flow

Configuration

Separate configuration avoids modifying pipeline source code

Orchestration

Nextflow runs the pipeline, submitting and monitoring jobs



Reproducible

Same results, every time




Portable

Write once, run anywhere



AT SCALE

Write once, scale everywhere



GROWTH

Nextflow usage continues to climb




Fig. 1 Google Scholar citation counts for bioinformatics workflow management systems. Sum of citations of the major publications of Galaxy, Nextflow, and Snakemake between 2018 and 2024 (Data in Additional File 1: Supp. Table 1)

DOI: 10.1186/s13059-025-03673-9



nf-core

The logo consists of the word "nf-core" in a bold, white, sans-serif font. The letter "x" in "core" is replaced by a stylized green and white apple icon with a stem and a single leaf.

A global community collaborating to build open-source Nextflow components and pipelines





A global community collaborating to build open-source Nextflow components and pipelines

140+

Pipelines

1,700+

Modules


13,600+

Community members

- Best-practice pipeline design patterns
- Standardized tooling, configuration & testing
- Extensive documentation for every pipeline



Community owned



Not just pipelines




Modules

Reusable process definitions for 1,700+ bioinformatics tools. Mix and match for custom workflows.



Tools

nf-core/tools CLI for creating, linting, testing, and launching pipelines with best practices.




```
nf-core pipelines create
```

The screenshot shows a terminal window with the title "nf-core pipelines create". The command "nf-core pipelines create" is typed into the terminal. The background of the slide features a teal vertical bar on the left side.

What's new in Nextflow

Releases 25.04, 25.10, and beyond



We've rewritten the heart of Nextflow

Nextflow is no longer a "Groovy DSL"

A new language parser means we can develop the language for users



Nextflow Language Server & VS Code Extension

A modern IDE experience for Nextflow development

Inline error reporting

Warnings and errors highlighted as you type

DAG preview

Visualize your workflow graph as you work

Quick navigation

Rename objects, click to jump straight to includes

Code completion

Context-aware suggestions at your cursor

Auto formatting

Format on save, with pre-commit or in CI

Compile time errors

Errors are detected at launch, rather than mid-run



DEVELOPER EXPERIENCE

Nextflow Language Server & VS Code Extension

A modern IDE experience for Nextflow development

```
workflow {
    def greetings = Channel.of('Ciao', 'Hello', 'Hola')
    def names = Channel.of('Dav `location` is not defined nextflow
    def locations = Channel.of(
        log.info("Launching workflow")
        sayHello(greetings, names, location)
        capitalize(sayHello.out) | view
}
```

[View Problem \(⌘F8\)](#) [Quick Fix... \(⌘.\)](#)



```
Channel.fromFilePairs(arg0) → Channel
```

Create a channel that emits all file pairs matching a glob pattern.

An optional closure can be used to customize the grouping strategy.

[Read more](#)

```
Channel.fromFilePairs(params.input)
```



workflows > methylseq > main.nf > workflow METHYLSEQ

```
// Aligner: bwameth
include { BWAMETH } from '../../../../../subworkflows/local/bwameth

/*
~~~~~ RUN MAIN WORKFLOW ~~~~~
*/
Preview DAG
WORKFLOW METHYLSEQ {
    take:
        samplesheet      // channel: samplesheet read in
        ch_versions      // channel: [ path(versions.yml)
        ch.fasta         // channel: path(genome.fasta)
        ch.fasta_index   // channel: path(star/index/)
        ch_bismark_index // channel: path(star/index/)
        ch_bwameth_index // channel: path(star/index/)


    main:
        ch_multiqc_files = Channel.empty()

        //
        // Branch channels from input samplesheet channel
        //
        samplesheet
            .branch { meta, fastqs ->
                single : fastqs.size() == 1
                return [ meta, fastqs.flatten() ]
            }

    emit
}
```

DAG Preview

The DAG diagram illustrates the data flow for the METHYLSEQ workflow. It starts with a 'samplesheet' channel, which branches into multiple paths. One path leads to 'CAT_FASTQ', which then connects to 'FASTQC'. Another path from 'samplesheet' goes to 'BISMARK'. From 'BISMARK', the flow splits into 'TRIMGALORE' and 'MULTIQC'. The 'MULTIQC' node also receives input from 'ch_versions'. Other nodes include 'QUALIMAP_BAMQC' and 'PRESEQ_LCEXTRAP', which also receive inputs from 'ch_versions'. The final output is 'multiqc_report'.



Static types & strict syntax

Nextflow's new native type system

Before

```
params.input = null

params.save_intermeds = false
params.threshold = 200
params.labels = 'one,two'
```

After

```
params {
    // No default value - required parameter
    input: Path

    save_intermeds: Boolean = false
    threshold: Integer = 200
    labels: String = 'one,two'
}
```



```
› nextflow run hello-error
```

```
N E X T F L O W ~ version 25.03.1-edge
```

```
Launching `hello-error/main.nf` [magical_euler] DSL2 - revision: 783f0990fe
```

```
ERROR ~ Script compilation error
```

```
- file : /Users/ewels/demo/nferror/hello-error/main.nf
- cause: Unexpected input: '{' @ line 15, column 10.
  workflow {
    ^
  }
```

```
1 error
```

```
NOTE: If this is the beginning of a process or workflow, there may be a syntax error in the body,  
such as a missing or extra comma, for which a more specific error message could not be produced.
```

```
-- Check '.nextflow.log' file for details
```



```
› NXF_SYNTAX_PARSER=v2 nextflow run hello-error
```

```
N E X T F L O W ~ version 25.04.1
```

```
Launching `hello-error/main.nf` [intergalactic_jepsen] DSL2 - revision: f22efb7ec8
```

```
Error main.nf:17:54: Unexpected input: '.'
```

```
| 17 |     Channel.of('Bonjour', 'Ciao', 'Hello', 'Hola') | .sayHello | view  
|
```

```
ERROR ~ Script compilation failed
```

```
-- Check '.nextflow.log' file for details
```



Check syntax ahead of time

```
> nextflow lint .
Error main.nf:21:1: Invalid include source: '/Users/ewels/GitHub/nextflow-io/rnaseq-nf/modules/dnaseq.nf'
| 21 | include { DNASEQ } from './modules/dnaseq'
| ^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^

Warn main.nf:37:3: Variable was declared but not used
| 37 | debug_variable = "Hello world"
| ^^^^^^^^^^^^^^

Nextflow linting complete!
✗ 1 file had 1 error
✓ 6 files had no errors
```



Workflow outputs & data lineage

Declarative output publishing and built-in provenance tracking

Workflow Outputs

Now out of preview in 25.10

- Declare outputs in a publish: section
- Assignment syntax for clear mapping
- Only in entry workflow — cleaner design
- Typed outputs with Channel annotations

Data Lineage

Built-in provenance (25.04+)

- Enable with: lineage.enabled = true
- Records every run, task, and output
- Explore lineage from the CLI
- Full audit trail for compliance





The complete platform for
scientific data analysis

Built by the creators of Nextflow



Launch, manage & monitor at scale



Launchpad

One-click pipeline launches across any infrastructure. Pre-configured parameters and compute environments.



Real-time monitoring

Live task tracking, resource metrics, cost estimation, and automatic failure notifications.



Governance & compliance

Complete audit trails, versioning, and role-based access control for regulated environments.




Multi-cloud

AWS, Google Cloud, Azure, and on-prem HPC. All from a single interface.

Node-RED + Seqera

Seqera Platform is your control plane for Nextflow



<https://seqeralabs.github.io/node-red-seqera/>

Cloud-native storage for Nextflow



Zero-copy I/O

POSIX interface to cloud storage. No "stage-in" copy.



Fusion Snapshots

Freeze task state on interruption, resume on new instance in < 2m



Multi-cloud

Works with all major cloud providers

SEQERA PLATFORM

And so much more..

The best place to run Nextflow

Seqera Studios

Interactive analysis environments on your compute

Wave: On-demand container provisioning

Builds optimized containers on the fly

Seqera Compute

Managed, optimized infrastructure

Seqera AI

AI-powered pipeline development

Seqera MCP, CLI, API, Terraform

Control your infrastructure however you want

Data Explorer

Browse your cloud data within Seqera



STARTED

Join the community

nextflow.io

Documentation, examples, resources

training.nextflow.io

Free, self-paced courses

community.seqera.io

Community forum

seqera.io/podcasts/

The Nextflow Podcast

seqera.io/blog/

The Open Science Blog

community.seqera.io

Community forum



Upcoming events



nf-core/hackathon

Online / distributed hackathon.

March 11-13



Nextflow Summit BOS

Boston Summit with talks,
training and hackathon. April 28



Nextflow Summit BCN

Online summit talks. Training +
hackathon in Barcelona. Oct '26

Plus many more: <https://seqera.io/events/upcoming/>



Thank you

Questions & Discussion

Phil Ewels

Senior Product Manager, Seqera

phil.ewels@seqera.io

nextflow.io · nf-co.re · seqera.io

